

Ca-pantothenate

Protocol number: M002236

Industry: Pharmaceuticals

Feed Size: < 300µm

Desired Fineness: d50 < 10µm

Quantity: ~ 200 g

Recommendation: Higher amount can be ground to d50 < 10µm with our Planetary Mill PULVERISETTE 5 classic line .

Result 1

Planetary Mill PULVERISETTE 5 classic line with 4 grinding bowl fasteners

main disk speed: 400 rpm

500 ml grinding bowl made of agate
+ 25x 20 mm agate grinding balls



Feed quantity: 200 g

Feed Size: <300 µm

Grinding time: 4 min

Final fineness: d90 < 63 µm

d50 < 15,5 µm

Comments: After 4 minutes of dry grinding, sample does not stick to bowl or balls. A longer dry grinding time is still possible (see following result). After 4 minutes of grinding, a fineness of d50 < 16 µm has been detected with our Laser Particle Sizer ANALYSETTE 22.

Result 2

Planetary Mill PULVERISETTE 5 classic line with 4 grinding bowl fasteners

main disk speed: 400 rpm

500 ml grinding bowl made of agate
+ 25x 20 mm Ø agate grinding balls



Feed quantity: 200 g

Feed Size: d₅₀ < 15,5 µm (see result 1)

Grinding time: + 2 min (Σ: 6 min)

Final fineness: d₉₀ < 59 µm
d₅₀ < 12,3 µm

Comments: After a total of 6 minutes of dry grinding, sample does not stick to bowl or balls. A longer dry grinding time is still possible (see following result). After 6 minutes of grinding, a fineness of d₅₀ < 13 µm has been detected with our Laser Particle Sizer ANALYSETTE 22.

A small portion of sample has been packed for demonstration.

Result 3

Planetary Mill PULVERISETTE 5 classic line with 4 grinding bowl fasteners

main disk speed: 400 rpm

500 ml grinding bowl made of agate
+ 25x 20 mm Ø agate grinding balls



Feed quantity: 180 g
Feed Size: d₅₀ < 12,3 µm (see result 2)
Grinding time: + 2 min (Σ: 8 min)
Final fineness: d₉₀ < 45,7 µm
d₅₀ < 9,8 µm

Comments:

After a total of 8 minutes of dry grinding, sample does not stick to bowl or balls. A longer dry grinding time is still possible (see following result). After 8 minutes of grinding, a fineness of d₅₀ < 10 µm has been detected with our Laser Particle Sizer ANALYSETTE 22.

Even longer grinding times are plausible to reduce the particle size.

Result 4

Planetary Mill PULVERISETTE 5 classic line with 4 grinding bowl fasteners

main disk speed: 400 rpm

500 ml grinding bowl made of agate
+ 25x 20 mm Ø agate grinding balls



Feed quantity: 180 g

Feed Size: d₅₀ < 9,8 µm (see result 3)

Grinding time: + 2 min (Σ: 10 min)

Final fineness: d₉₀ < 36,2 µm
d₅₀ < 8,7 µm

Comments: After a total dry grinding time of 10 minutes, sample starts sticking lightly to bowl and balls.
For a further reduction of particle size, we recommend grinding in suspension.

The present d₅₀ < 8,7 µm has been detected with our Laser Particle Sizer ANALYSETTE 22.

Contact to our Applications Laboratory: Leos Benes · Phone: 0049 67 84 70 122 · benes@fritsch.de

